



Specification

1. Power Source : One 1.5V AA Battery
2. Laser Output : 6mW (LLLT, Low Level Laser Therapy)
3. Laser Wavelength : 635nm
4. Net Size : Height 210mm (8.27 inch) * Diameter 24mm (0.95 inch)
5. Case Size : 242mm(H) * 85mm(W) * 50mm(D) (9.53 inch * 3.35 inch * 1.97 inch)
6. Net Weight : 80g(2.82oz) (Battery included)
7. Gross Weight : 280g(9.88oz)
8. Type of bristle : Ultra soft and thin bristles
9. Color Variation : Silver, Gold, Green, Violet and Blue (All Metallic)
10. Origin : Made in Korea
11. Warranty : 1-year warranty against product defects.
12. Certification : CE 1023, ISO 13485, CB, CE, SASO, GOST, RoHS, KFDA, GMP



M&H (Medical & Human Technology)
 436-16, Changgok-ri, Paltan-myeon, Hwaseong-si, Gyeonggi-do, Korea
 Tel : +82-31-452-1495 Fax : +82-31-452-1497
 david@mnhtw.com http://www.mnhtw.com



Component

1. Laser Toothbrush
2. 2 replacement brush heads (large and small)
3. One 1.5V AA Battery
4. Case
5. Manual



What is Dr. M[®]?

It is the world's first laser toothbrush that uses low level laser therapy to treat dentin hypersensitivity (sensitive teeth). Dr. M[®] is more than just a pain reliever, it cures for sensitive teeth. It is clinically proven by leading hospitals and universities. It's a medical device.

What causes sensitivity?

As we age several things happen to our teeth. First, a lot of people grind their teeth so they are losing enamel on teeth. However, the biggest one is the recession of gums it can happen from brushing too hard. Do not forget always to use soft toothbrush. It can be as bad as not being able to eat ice cream, not being able to drink cold beverages or even hot water and not even going to get teeth clean because you are so sensitive when you do that.

How does Dr. M[®] cure the hypersensitivity?

Dentin hypersensitivity in most cases is related to exposure of root surface by recession of gums. A tooth with dentin hypersensitivity has dentinal tubules twice the size of a normal tooth's dentinal tubules. The purpose of Dr. M[®] is to treat hypersensitivity **by reducing the enlarged diameter of the affected dentinal tubules.**

How long does it take for the treatment?

Each treatment takes five minutes while brushing your teeth. For the best results, it is recommended to use Dr. M[®] on problem area two to three times per day over a period of four to six weeks.

Operation



M&H (Medical & Human Technology)
 436-16, Changgok-ri, Paltan-myeon, Hwaseong-si, Gyeonggi-do, Korea
 Tel : +82-31-452-1495 Fax : +82-31-452-1497
 david@mnhtw.com http://www.mnhtw.com



Specialty

1. Laser is only emitting when the brush is in the mouth.
2. Biological safety test on bristles and brush head : Cytotoxicity (ISO 10993-5), Sensitizing Test (ISO 10993-10) and Intracutaneous(intradermal) reactivity (ISO 10993-10)
3. Ultra soft and thin bristles
4. After 5-minute operation it automatically shuts off its power.
5. Nano-Silver is coated on the surface where bristles are located. (Anti-bacterial function)
6. Sensor is coated with gold.
7. Ultra thin and soft bristles help you to have cleaner teeth.
8. It really works, clinically proven. Dentists recommended.

